

E CATALOG

---

**PLASMA STERILIZATION  
DEVICE FOR GENERAL  
PURPOSE**

---

**HRF 3000**

**Hydrogen Peroxide  
(H<sub>2</sub>O<sub>2</sub>) Plasma  
Sterilization Device**

---

**PRODUCER:** TEKNOMAR MAKİNA İMALAT  
İTHALAT İHRACAT SANAYİ ve TİCARET LTD.  
ŞTİ.

**Notified Body No:** CE-1984

---

**THE MOST EFFECTIVE AND SAFE SOLUTION  
IN STERILIZATION**

4th generation HRF COLD PLASMA Technology

**"Sterilization will never be interrupted"**

## HRF 3000 H<sub>2</sub>O<sub>2</sub> Plasma Device Automatically Detects Damp Load, Dries and Sterilization Process Continues.

### High Technological Feature:

#### RELIABLE AND ECOLOGICAL PERIOD IN STERILIZATION WITH HRF 3000

HRF 3000 Hydrogen Peroxide Sterilization Device; All kinds of plastic, polymer, inorganic, electromechanical instruments and surgical instruments, etc .. re-used, all kinds of heat and moisture sensitive medical equipment that sterilizes the new generation of ecological hydrogen peroxide plasma sterilization device.

HRF 3000 (H<sub>2</sub>O<sub>2</sub>) Hydrogen Peroxide Plasma Sterilization Device is a low temperature, moisture free, non-toxic sterilization device.

Used for sterilization of all kinds of medical equipment and medical instruments such as metal and non-metal, inorganic, polymer, medical grade, etc.

It has **HRF Cold Plasma** system which converts hydrogen peroxide used during sterilization into water and oxygen.



It does not require any installation and only electrical connection. Sterilization temperature can be set at 35 - 55 and optionally at intermediate values. Works with 60% H<sub>2</sub>O<sub>2</sub> cartridge.

After sterilization, the materials are presented ready to use without the need for extra ventilation with device's technological design and applications.

No harmful waste, the final products are water vapor and oxygen. Therefore, it provides safe sterilization as well as safe use for personnel and the environment.

## HRF 3000 Hydrogen Peroxide Gas Plasma Technology Advantages and Disadvantages

- Extends the life of medical devices and instruments.
- Cost effective
- Shortens the processing time in the sterilization center (CSSD)
- Ecological.
- 37 - 55 C Operating Temperature
- Plasma in chamber
- Tyvek® package compatible.
- No toxic waste.
- Does not need any ventilation line
- User and environment friendly
- 10<sup>-6</sup> (Sterility Assurance Level)
- Sterilization is not addicted to packaging consumables.



### Disadvantages


- Materials that absorb hydrogen peroxide, cause catalytic degradation of hydrogen peroxide, react with hydrogen peroxide materials such as organic sulfides, silver are not suitable for sterilization.
- No suitable sterilization method for cellulose, fabrics and liquids,
- Use of cellulose-free polypropylene, synthetic material such as Tyvec for packaging
- **Controlled use is recommended for sterilization of lumen materials.**

### Packaging Specifications:

Cellulose-free, synthetic materials such as polypropylene or Tyvec

## HRF 3000 STERILIZATION PROCESS

- **Preconditioning** (Vacuum Test, conditioning, ... etc.),
- **Sterilization** (H<sub>2</sub>O<sub>2</sub> injection, sterilization, diffusion, plasma, ... etc.),
- **Aeration** (Air Washing),
- All stages start and end automatically without user intervention.
- If there is H<sub>2</sub>O<sub>2</sub> in the chamber, the doors will not be opened for safety reasons without air washing.



HRF 3000 Lumen Full Cycle Sterilization Diagram


## HRF 3000

In sterilization safety, the sterilization device uses a fixed programmed sterilization cycle that eliminates potential user error.

Provides simple and fast operation with predetermined program options for language selection and sterilization on the large touch screen.

**Complies with ISO 14937 Sterilization validation standards.**

There are 3 H<sub>2</sub>O<sub>2</sub> Sterilization Programs including Short, Long and Lumen.



**Lumen Sterilization / (PCD)**

HRF 3000 (H<sub>2</sub>O<sub>2</sub>) Hydrogen Peroxide Plasma Sterilization Device performs the sterilization process of the following types of endoscope, device, instruments, disposable materials, medical products and other similar materials.

Does not damage plastic and electronic equipment due to low temperature sterilization, instruments such as camera head, fiber optic cable, rigid and flexible endoscopes are safely sterilized.

| <b>ENDOSCOPES</b>  |  |
|--|--|
| <b>Controlled use is recommended for sterilization of Multi-Complex lumen materials.</b>   |  |
| <b>Rigid Endoscopes / lumens</b>   | <b>Flexible Endoscopes / lumens</b>  |
| Laryngoscope<br>Arthroscopy<br>Laparoscopes<br>Trocar Cannula<br>Trocar Case<br>Resectoscope<br>And such   | Bronchoscopy<br>Ureteroscopy<br>Hysteroscopy<br>Cystoscopy<br>Koledoskop<br>And such   |
| <b>Device / Instrument</b>   |  |
| implants<br>Defibrillator Pedals<br>Electrocautery Products<br>Esophageal Dilators<br>Kri-Probe<br>Doppler<br>Head Pressure Transducer Cables<br>Endoscopic Products and etc | Fiber Optic Cables<br>Laser Hand Products<br>Fiber Accessories<br>Ophthalmic Lenses<br>Radiation Therapy Equipment<br>Surgical Power Equipment<br>Drilling Tools<br>Ultrasound Probes<br>Video Camera and connection apparatus and etc |



**FTPE 850mm ø1mm**

**FTPE 1000mm ø2mm**

**FTPE 1200mm ø2mm**

**Teknik Özellikler / Kısa Datasheet**

| MODEL   | S75                 | S125                | S175                |
|---|---------------------|---------------------|---------------------|
| TYPE  | HRF 3000            | HRF 3000            | HRF 3000            |
| SCREEN  | Touch screen 7"     | Touch screen 7"     | Touch screen 7"     |
| OUTER DIMENSIONS (W-D-H) (mm)                             | 730x955x1860        | 730x955x1860        | 820x955x1860        |
| INNER DIMENSIONS (W-D-H) (mm)                             | 440x370x700         | 440x500x700         | 550x500x700         |
| CHAMBER VOLUME (lt)                                       | 113                 | 154                 | 192                 |
| EFFECTIVE VOLUME OF CHAMBER (lt)                          | 92                  | 131                 | 171                 |
| CARTRIDGE H <sub>2</sub> O <sub>2</sub> SOLUTION QUANTITY | 120 cc              | 120 cc              | 120 cc              |
| *POWER (WATT)(Max.)                                       | 3200W               | 3200W               | 3200W               |
| WORKING VOLTAGE   | 230 VAC, 50/60 Hz   | 230 VAC, 50/60 Hz   | 230 VAC, 50/60 Hz   |
| EXTERNAL CHASSIS MATERIAL                                 | STAINLESS STEEL 304 | STAINLESS STEEL 304 | STAINLESS STEEL 304 |
| CHAMBER MATERIAL  | STAINLESS STEEL 316 | STAINLESS STEEL 316 | STAINLESS STEEL 316 |

All models have double door selection (D75,D125,D175)  
For more information check the Datasheet.

| HRF 3000   | CERTIFICATES            |
|--|-------------------------|
| CE Certificate- 93/42/EEC                          | YES                     |
| EN ISO 13485                                       | YES                     |
| EN ISO 9001  | YES                     |
| TUR (Technological Product Experience Certificate) | YES                     |
| Domestic Goods Certificate                         | YES                     |
| Free Sales Certificate                             | YES                     |
| UTS Registered                                     | YES                     |
| TEYDEP Project Success Certificate                 | YES                     |
| LVD / EMC Tests and Certificates                   | YES / EN 60601-1-1:2009 |
| EN ISO 14937 Certificates                          | YES                     |
| Prion Tests  | YES                     |
| Lumen Material Sterilization Test                  | YES                     |
| Residue Test                                       | YES                     |
| Corrosion Test                                     | YES                     |
| Device Type Test                                   | YES                     |

| Ethylene Oxide And Hydrogen Peroxide Sterilization Comparison Table |                                      |  |
|---|--------------------------------------|--|
| Sterilization Method  | Ethylene Oxide                       | Hydrogen Peroxide                                      |
| Teknomar  | EO - C <sub>2</sub> H <sub>4</sub> O | H <sub>2</sub> O <sub>2</sub>                          |
| Sterilization Time  | ~ 4-12 Hours                         | ~ 70 Minutes   |
| Sterilization Cost  | ~ 12, <sup>00</sup> * €              | ~9, <sup>00</sup> * €                                  |
| Cartridge Storage Condition   | Compelling                           | Appropriate  |
| Danger Level  | High                                 | Ecofriendly  |
| Level of Preparation for Sterilization                              | Medium                               | Long   |
| Capacity  | High Capacity                        | Limited Capacity                                       |
| Material Compatibility  | Mixed Material                       | Different Programs for Different Products              |
| Air and Environment Pollution                                       | Risky                                | Ecological   |
| External Connection   | Necessary                            | Unnecessary  |
| Operating temperature   | 37-55 °C                             | 37-55 °C   |
| Sterility Assurance (SAL)   | 10 <sup>-6</sup>                     | 10 <sup>-6</sup>                                       |
| Maintenance   | Expensive                            | Reasonable   |
| Installation  | Hard                                 | Easy   |
| Working Principle   | Only EO                              | H <sub>2</sub> O <sub>2</sub> Gas Vapor or Plasma Only |
| Diameter and Length   | Unlimited                            | 1 mm Ø - 850 mm / 2 mm Ø 1200 mm                       |
| Residue on Product  | Risky                                | No residue on product / partial lumen                  |
| Approximate Cost of Device  | ~ 20.000,00 €*                       | ~ 40.000,00 €*   |
| Total Score   | <b>6/17</b>                          | <b>8/17</b>  |

The table shows the comparison of sterilizers in terms of different parameters and properties. Although sterilization is precise in each of the methods, there are variables that must be considered from the end-user point of view. The ranking shows that the most suitable sterilizer is the 8/17 rated Low Temperature & Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) Gas Plasma Sterilization Device

E CATALOG

PLASMA STERILIZATION  
DEVICE FOR GENERAL  
PURPOSE

HRF 3000

Hydrogen Peroxide  
(H<sub>2</sub>O<sub>2</sub>) Plasma  
Sterilization Device

**TEKNOMAR®**

**PRODUCER:** TEKNOMAR MAKİNA İMALAT  
İTHALAT İHRACAT SANAYİ ve TİCARET LTD.  
ŞTİ.

**Notified Body No:** CE-1984

TEKNOMAR MAKİNA İMALAT İTH. İHR. SAN. VE TİC. LTD.ŞTİ.

Ostim OSB Mahallesi 1269. Cadde No: 29

Yenimahalle ANKARA / TÜRKİYE 06370

Tel: +90 312 385 00 40

Fax: +90 312 385 67 84

[www.teknomar.com.tr](http://www.teknomar.com.tr)

[info@teknomar.com.tr](mailto:info@teknomar.com.tr)